


The Simple Solution to Zeta Potential Measurements  
in Diluted Colloids


# Zeta Potential Meter

Model 502



## Rotating Prism Technique

When light passes a rectangular prism, the object seems to move as a course of the light is out of position when a prism was rotated slowly. A microscope has a prism (rotating prism) in conjunction with a galvanometer built-in, and the rotating speed and rotating direction of prism can be regulated. The particle is observed with a television monitor with a lattice. Many particles move electrophoretically when we apply the electric field. The zeta potential value of the particle (20 degrees Celsius, water system) is displayed digitally when we adjusted the observed particles to appear stationary using the rotating prism technique.



**Nihon Rufuto Co., Ltd.**

## FEATURE

### •Traditional measurement principle

Model 502 measures the Zeta potential of colloidal particles by traditional MicroElectrophoresis. There is no black box compared with the laser scattering electrophoresis.

### •Observation of nanoparticles by laser dark field illumination

502's unique ultra-microscope design provides a high contrast image even for nano-sized particles. Depending on the refractive index of the particles relative to the suspending media, measurements can be made on particles as small as 20nm.



Picture of 130 nm PSL

### •Observation of the particles in stationary layer

High level laser optical system and high performance CCD camera can observe only the particles in stationary layer. It provides a highly precise result without being affected by the electroosmotic flow.

### •The rectangular electrophoresis cell which was put up horizontally

The horizontally stood rectangular cell provides a highly precise result. Because the layout is hard to occur the non-symmetrization of the electroosmotic flow caused by the sedimentation of particles on the chamber wall.

## TECHNICAL DATA

<p><b>Principle:</b> User adjust rotating prism for stationary image. Or measurement of the transit time of particle by stopwatch</p> <p><b>Measurement range:</b> <math>\pm 100</math> mV</p> <p><b>Particle size range:</b> Typically 40 nm to 50 <math>\mu</math> m depending on refractive index and sedimentation</p> <p><b>Sample required:</b> 8 cc</p> <p><b>Optical system:</b>            Light source: 632.8 nm He-Ne laser            Microscope magnification: <math>\times 280</math>            TV monitor: CCD monochrome camera and 8.4" monochrome LCD monitor</p>	<p><b>Electrophoresis chamber:</b>            cross section: 1 <math>\times</math> 10 mm            construction: pure silica            electrode separation: 4.88 cm            max field strength: 30 V/cm            anode: molybdenum            cathode: palladium            electrode compartments: derlin</p> <p><b>Supply voltage:</b> 0 ~ 150 V</p> <p><b>Power required:</b> 100 VAC, 50/60 Hz</p> <p><b>Dimensions:</b></p> <p><b>Weight:</b></p>
--	--

*Specifications and descriptions in this brochure subject to change without notice.*



**Nihon Rufuto Co., LTD.**  
 Scientific Instruments Division

Uenofuji Bldg. 5-1-8 Higashiueno  
 Taito-ku, Tokyo 110-0015 JAPAN  
 TEL:+81-(0)3-3847-6880  
 FAX:+81-(0)3-3847-6890  
<http://www.nihon-rufuto.com>

Distributor

KS110418