

The logo for Occhio, featuring the word "Occhio" in a white serif font, with a stylized eye graphic to its right. The background is a dark blue gradient.

Occhio

You need precision, you want morphology



OCCHIO Flow-Cell FC200S+

The best solution for measuring suspensions, emulsions and foams





Imaging solutions in particle analysis

www.occhio.be



OCCHIO Flow-Cell FC200S+

Size , shape, counting for suspensions, emulsions and foams



Through the efforts of an international and multidisciplinary team of engineers, **OCCHIO** offers you a complete range of solutions, starting from 200 nanometers and ranging up to centimeters.

Whether it is for laboratory instrumentation, «at line» or even «on line» solutions, **OCCHIO** is prepared to be your partner in high-level powder characterization. **OCCHIO** and **OCCHIO Flow-Cell** bring you accuracy, profit and innovation.

_ Accuracy


With its proprietary Light and high quality lens, **OCCHIO Flow-Cell** will change your own perception of image analysis, measuring suspensions or emulsions which are invisible under normal microscopy.

_ Profit

OCCHIO Flow-Cell is an automatic device dedicated to emulsions or suspensions quality characterization. It is easy to use and carries out rapid analyses in less than few minutes.

_ Innovation

Morphology measurement is more than shape description. To improve, you need robust and significant measurement. Based on decades of university research, the **OCCHIO Flow-Cell** provides your R&D and Production departments with dedicated parameters, specially engineered for your industrial purposes.



_ Size measurements (from 200 nm up to 1000 µm)

Area diameter - Mean diameter - Length - Width - Maximum distance

_ Shape parameters

Elongation - Circularity - Convexity - Shape factor - Luminance & Special parameters.....

For more information contact us

info@occhio.be www.occhio.be